




exploded axon of interior medical capsule

exploded axon of BlueCore insert showing the cellular divisions of each segment.


Interior capsule is a soft insert made of ripstop nylon. Sections zip away from each other for flexible distribution and use of medical kits.

BlueCore medical insert is designed to be pulled out of the pack quickly in the case of an emergency. The user flips the lid open and pulls insert up and out of the pack.



base segment doubles as wash basin


ripstop nylon segments zip together to form a compressible cylindrical tower



insert fits into the core of the pack interior compression straps secure insert to aluminum stays

BlueCore is packed with the lightest gear towards the outside of the pack and the heaviest gear towards the frame. Pack lid opens to top-loading body.

When the lid of the pack is removed, the pack's components are easy to reach from the top entry.



lid holds quick access items such as knife, snacks, map, stakes and headlamp

BlueCore interior medical capsule

clothing and personal items

cookware, food, water pump

sleeping bag, tent with rainfly and groundmat in compression sacks

personal hygiene items

mat, tent poles




4" ace wraps




bandages



exam gloves



6 oz. medicine bottles



backcountry gear (listed above)

BlueCore interior medical capsule

lateral storage zones for water and carryout

ergonomic back support

hybrid interior-exterior perimeter frame

